



WALL MOUNTED UNIT WITH SPECULUM BOX

WALL MOUNTED UNIT



UNIT

- Operates with 220V 50/60Hz mains power.
- Comprises two handles with adjustable light intensity.
- The handles can extend up to 3 meters thanks to the spiral, durable cable.
- The handle automatically turns on when removed and turns off when replaced by its unique magnetic switching technology.
- The automatic turn-off function turns off the light on the headrests when the handle is left on for more than 5 minutes.
- Fast soft start and soft end function extend the lamp's life.
- Stylish, illuminated, and waterproof on/off switch.
- Connector lock system used for easy instrument installation.
- Separate units can be added if required



SPECULUM BOX

- *It has 4 compartments.
- *Speculas of various diameters can be used.
- *The mounting brackets allow easy and convenient installation.
- *Each compartment holds 15 speculums, for a total of 60.